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Testo

DESTINAZIONE D'USO

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USER INSTRUCTIONS
Text

INTENDED USE

Il Kit **NtDrain** is a kit used to clamp a **ventricular drain** catheter comprising a device implanted in the hole drilled in the cranial bone and specific disposable instruments. The device prevents unwanted movement of the catheter from its position. After the catheter is removed, sliding it out as a closed procedure, the **NtDrain** device will remain permanently implanted in the cranial vault to close the drilled hole, thereby avoiding the formation of a hollow in the skin with evident aesthetic benefits. The **NtDrain** device is used to **lock catheters having external diameters of 3.3 and 3.0 mm**.


INTRODUCTION

- The **NtDrain** device must not be implanted if patients are under the age of 16 years.
- The **NtDrain** device should only be implanted by neurosurgeons who are fully aware of the intended purpose, surgical technique, the specific instruments and the risks arising from improper use.
- The device penetrates the skull by 3.5 mm. The neurosurgeon verifies, in cases of particularly thin cranial theca, that **NtDrain** maintains the necessary distances from the surface of the dura mater.
- The **NtDrain** device must be implanted in an operating room under aseptic conditions.

PRODUCT INFORMATION


The **NtDrain** kit comprises:

- **NtDrain** made up of:
 - a **crown** (fig.1) with external teeth designed to grip the craniotomy hole and internal flexible tabs used to clamp the catheter
 - a central, internal **insert** (fig.2) designed to widen the flexible tabs so that the catheter can be released.
- A disposable polymer applicator **NtApplier** (fig.3) for the NtDrain implant
- A disposable polymer activator **NtActivator** (fig.4) for extraction, from under the skin, of the NtDrain catheter.

The device is inserted in the hole drilled into the cranial vault, it can block/release the catheter during the insertion stage using **NtApplier**. When no longer needed, the catheter can be removed using **NtActivator** acting externally on the skin by pressing the **insert** positioned in the **crown**. The **NtDrain** kit is packaged in 2 rigid, sterile double blisters in a cardboard box with two independent zones, containing:

- NtDrain pre-assembled on NtApplier (zona 1)
- NtActivator (zona 2)
- User instructions
- Adhesive traceability labels. (Zona 1)

WARNING: NtDrain, is made with radiolucent Peek OPTIMA and does not interfere with CT and MRI examinations.

WARNING: This is a **single-use** disposable device. Do not sterilise again and/or reuse. Reuse of the device may entail the risk of malfunction, breakage, or infection. The device is suitable for long-term insertion.

WARNING: Do not use the **NtDrain** device with patients who have shown signs of allergic reaction to PEEK.

WARNING: The **NtDrain** device is available as two variants specifically for **3.3 and 3.0 mm catheters**. The hole drilled in the skull must be made using a **14 mm diameter drill**.

INFORMATION ABOUT INSTRUMENTATION

The **NtDrain** device is equipped with the following dedicated instruments:

- **APPLIER FORCEPS** (NtApplier) fig.3
- **ACTIVATOR** (NtActivator) fig.4

NtApplier is pre-fitted with NtDrain and is used to slide the device along the catheter and facilitate insertion into the craniotomy hole. The handle of the forceps has dedicated geometry to open, by means of compression on the insert, the tabs of the crown.

NtActivator is used to release the device from the catheter by opening the internal tabs of the crown by compressing the insert under the skin.

SURGICAL TECHNIQUE

NtDrain device implantation

The **NtDrain** device must be implanted in a craniotomy hole made using a **14 mm drill**. Decisions concerning insertion technique, verification of correct positioning and the choice of catheter type for intracranial drainage must be taken by the surgeon on a case by case basis. The selection of eligible patients is the responsibility of the surgeon.

WARNING: Do not use the **NtDrain** device with patients who have shown signs of allergic reaction to PEEK.

WARNING: The **NtDrain** device is available as two variants specifically for **3.3 and 3.0 mm catheters**. The hole drilled in the skull must be made using a **14 mm diameter drill**.


ATTENTION:

- Make sure the package is in good condition. In the case of damaged packaging or packaging with labels that are no longer legible, do not use the contents.
- The **NtDrain** device must be used before the expiry date indicated on the packaging.
- After opening the 2 secondary blisters, make sure that the **2 NtDrain** device is not defective or damaged, in which case it must not be used: the manufacturer must be notified.

Insertion stage:

The surgeon

- must prepare a skin flap suitable for cranial perforation using a 14 mm drill, making sure that the incision line does not coincide/does overlap the hole even partially.

code	mm
0202-043.0	3,0
0202-043.3	3,3


CATHETER PLACEMENT

- (zone 1) Using **NtApplier**, already pre-assembled with **NtDrain**, the surgeon can compress the insert and open the tabs of the crown to facilitate insertion of the catheter, with its stiletto fitted, into the **NtDrain** insert.
- Once the catheter has been inserted with its stiletto reaching the ventricle at the required depth, the surgeon can secure it with ordinary surgical forceps and then slide **NtDrain** along the catheter to compress the insert and open the tabs.
- The surgeon will then be able to slide it so that the crown of the device is inserted into the hole drilled in the bone, remove the dedicated forceps and, with finger pressure on **NtDrain**, insert it into the cranial hole. The device will stabilize thanks to the teeth on its crown.
- The **handle** of **NtApplier** can also be used to vary the depth of the catheter by pressing more or less on the insert of the device that opens or closes the flexible tabs that release or tighten the catheter.
- Lastly, the catheter is locked by the flexible tabs and it is possible to remove the stiletto and proceed with the tunnelling under the skin and suturing the skin flap.

LOCKING MECHANISM

ATTENTION: Before carrying out the tunnelling of the catheter, check that the catheter itself is positioned in the recess on the NtDrain insert.

(Step 1)



NtDrain can be removed in the event of **difficulty**, by removing the catheter by compressing the device's **insert** using its **handle** to open the **flexible tabs**. Remove the device using a dissector and by levering on the crown.

Catheter removal stage: (zone 2)

The surgeon

- must palpate the skin flap to identify the position of the catheter under the skin.
- Using **NtActivator** positioned on the swelling indicating the catheter emerging from the hole, compress the skin in a direction at right angles to the underlying **insert**; this will allow the **flexible** **tabs** to open and release the catheter (**Step 2**).
- Maintaining compression on the insert to keep the **flexible tabs** open, the catheter can be easily removed (**Step 2**).

CATHETER REMOVAL

NtDrain remains under the skin and, after fibro-integration, will perform its **aesthetic function**.

Kit per fissaggio catetere con diametro esterno 3,0 mm composto da NtDrain 3,0, NtApplier, NtActivator	Kit for fastening a catheter with an external diameter of 3,0 mm comprising NtDrain 3,0, NtApplier, NtActivator
Kit de fixation d'un diamètre externe de 3,0 mm, composé par NtDrain 3,0, NtApplier, NtActivator	Kit de fixation d'un diamètre externe de 3,0 mm, composé par NtDrain 3,0, NtApplier, NtActivator
Kit para fijación del catéter con diámetro exterior de 3,0 mm, compuesto por NtDrain 3,0, NtApplier, NtActivator	Kit para fijación del catéter con diámetro exterior de 3,0 mm, compuesto por NtDrain 3,0, NtApplier, NtActivator
для фиксации катетера наружным диаметром 3,0 мм, состоящий из NtDrain 3,0, NtApplier, NtActivator	для фиксации катетера наружным диаметром 3,0 мм, состоящий из NtDrain 3,0, NtApplier, NtActivator

KIT NtDrain 3.0 Code 0202-043.0 KIT NtDrain 3.3 Code 0202-043.3

DISPOSITIVO DI BLOCCAGGIO CATETERE – ISTRUZIONE PER L'USO
DISPOSITIF DE BLOCAGE DE CATÉTER – MODE D'EMPLOI
DISPOSITIVO DE BLOQUEO DEL CATÉTER – INSTRUCCIONES DE USO
УСТРОЙСТВО ДЛЯ ЗАКЛИМА КАТЕТЕРА – ИНСТРУКЦИЯ ПО ЭКСПЛУАТАЦИИ

DISPOSITIVO DI BLOCCAGGIO CATETERE – ISTRUZIONE PER L'USO
DISPOSITIF DE BLOCAGE DE CATÉTER – MODE D'EMPLOI
DISPOSITIVO DE BLOCAGE DEL CATÉTER – INSTRUCCIONES DE USO
УСТРОЙСТВО ДЛЯ ЗАКЛИМА КАТЕТЕРА – ИНСТРУКЦИЯ ПО ЭКСПЛУАТАЦИИ

ntplast®
advanced biomedical solutions

NTPLAST s.r.l.
Via G. Donizetti, 20 00198
Roma - ITALY
info@ntplast.com - www.ntplast.com