ntplast

ntcranial small

Istruzioni per l’uso_ rev. 03 / 2026

User instructions_ rev. 03 / 2026

Bedienungsanleitung_ rev. 03 / 2026

Instrucciones de uso_ rev. 03 / 2026

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ENGLISH

Intended use

ntcranial small is a plastic device intended for the stable fixing of the craniotomy flap (variant of **ntcranial**) to be positioned only along the craniotomy cut. It consists of a tensor and a plate.

Suitable for both autologous bone and synthetic cranioplasties.

Introduction

- The **ntcranial** small device should only be implanted by neurosurgeons who are fully aware of the intended use, surgical technique, risks arising from improper use, and the instruments dedicated to the device.
- The **ntcranial** small device must be implanted in an operating room under aseptic conditions.
- When implanting **ntcranial** small devices, the utmost attention must be paid to the movements.
- The use of **ntcranial** small devices combined with other cranial fixing devices of competitors is not recommended. If this is the case, it is the surgeon's responsibility to verify the stability of the fixed flap.

Product information

The **ntcranial** small (variant of **ntcranial**) is a long-term implantable device.
The **ntcranial** small device consists of a tensor a locking plate (fig. 1).
The lower part of the tensor will have bone elastic arms with a new design. The transverse arms are no longer cylindrical but oval, slightly thinner, with rounded ends and a shape that helps slide them into the epidural space.
The cylindrical plate, assembled on the upper part of the tensor, will be equipped with elastic hooks that can descend along the teeth of the tensor.
The shape of the device allows for total contact of the plate with the transverse arms thus allowing tightening of the device on bones of minimum size.
After assessing the shape and size of the flap, a surgeon can decide the number of devices and choose between the two formats, i.e., **ntcranial** and **ntcranial small**.
ATTENTION: The **ntcranial** small device is indicated for bone flaps ranging between 0.5 mm and 20 mm in thickness.
ntcranial small is suitable for positioning in craniotomy cuts with a width of between 1.7 mm and that resulting from the use of standard craniotome cutters.
If ntcranial small devices are used in conditions other than those indicated above, fixing of the flap may prove to be inadequate, as it may be unstable and misaligned with respect to the cranial vault.
MATERIAL: The **ntcranial** small device is made entirely of biocompatible PEEK OPTIMA® plastic, which does not interfere with CT or MRI exams.
WARNING: Do not use the **ntcranial** small device on patients who have shown signs of allergic reaction to PEEK.

Sterilisation

The device (Fig. 2) is packaged in double blister packs and supplied sterile.
WARNING: The ntcranial small device is for single-use only. Do not sterilise again and/or reuse.
Reuse of the device may entail the risk of malfunction, breakage, or infection.

Package Contents

The packaging of the **ntcranial** small device consists of an outer cardboard box containing 1 sterile device, already assembled in double semi-rigid medical paper blister with 321 mm FITFORM 200 F film. The cardboard box also contains the use **instructions** and **3 labels for traceability**, attached with double-sided adhesive, to each double blister.

Contraindications

Contraindications include, but are not limited to, the following:

- hypersensitivity or allergies to PEEK
- latent or active intracranial infections
- fever and leukocytosis
- bone conditions that are not suitable for the application of **ntcranial small** device
- bone infections
- subtentorial craniotomy
- post-traumatic multi-fragmented bone flaps
- patients under 10 years.

Complications

Improper use may cause the breakage of the **ntcranial** small device or instability during fixing, with possible dural or brain damage.
The general perioperative risks include:

- infections,
- brain injury, neurological deficit
- haematomas at the implantation site
- alterations in wound healing.

Surgical technique

The **ntcranial** small device must be implanted following only the instructions in these instructions.
The identification of the number of ntcranial devices and their configurations together with the use of the ntcranial small falls under the sole responsibility of the surgeon, as does the choice of the location and the correct application of the devices.
The selection of eligible patients is the responsibility of the surgeon.

Positioning the devices

- Prepare, at the surgeon's request, the number of packaged devices needed to fix the flap.
- Make sure the packages are intact.** In the case of damaged packaging or packaging with labels that are no longer legible, do not use the contents. The device must be implanted before the expiration date indicated on the packaging. Check also that the **ntcranial Cutter** is intact.
- Remove the device from the packaging without compromising its sterility. **If any component appears to be defective or damaged, do not use the device.**
- ntcranial small** is already assembled (fig. 3).
- WARNING: Do not try to alter the curve of the elastic arms of the tensor by acting directly on it with your hands or tools (e.g., pliers).** Altering the elastic arm curve may hinder proper assembly with the bases. **WARNING:** Place **ntcranial small** along the craniotomy cut in the extradural, subdtecal site (Fig. 4).
- WARNING:** Make sure that the **ntcranial small** is correctly oriented (Fig. 5).
- WARNING: Do not place the ntcranial small at the perforator hole.** If the device is incorrectly positioned, locking the device may be impossible or unstable.
- It is recommended that ntcranial and ntcranial small devices should be positioned as evenly as possible along the perimeter of the craniotomy: do not neglect any side of the craniotomy** (Fig. 4). The number and type of devices must be compatible with the shape and area of the flap. **Do not use less than 3 devices to achieve stable fixing.** For standard craniotomies, 3 or 4 devices are generally sufficient. More devices may be required for large craniotomies.
- After inserting all the **ntcranial** and **ntcranial small** devices, place the bone flap in the craniotomy site.
- While holding the tensor in place with two fingers, bring the locking plate into contact with the bone by pushing it down with the fingers of your other hand (Fig. 6). Repeat this step for all the devices (Fig. 7).
- WARNING: When lowering the plates, it is important to make sure that the bone flap is kept in the correct position.**
- Take the tensor again with two fingers, pulling it upwards to activate the underlying elastic component. At the same time, push the plate down as far as possible using two fingers of your other hand. Compared to the previous step, the locking plate will descend further, consolidating the fixing. Repeat this step for all the **ntcranial** and **ntcranial small** devices.
- ATTENTION: Perform the steps above by hand only. It is absolutely forbidden to use tools (such as pliers) that could damage the ntcranial and ntcranial small devices.**
- Once all the **ntcranial** and **ntcranial** small devices are positioned, check the fit and stability of the flap.
- 10) Remove excess tensor lengths by means of the supplied **ntcranial Cutter** (Fig. 8). The cutting of the tensor must be flush with the upper surface of the locking plate.

ATTENTION: During cutting, hold the handle of the tensor with your fingers so that it does not come off accidentally (Fig. 8).

WARNING: To remove the tensor parts, it is recommended to use only the **ntcranial Cutter**.

WARNING: Check that the resulting cut profile is smooth. If necessary, trim. Be careful when doing this: small burrs may pierce the surgeon's gloves.

Removing the devices

Using the **ntcranial Cutter**, cut the locking plate (Fig. 10). Repeat this step for each **ntcranial** and **ntcranial** small devices.

Post-operative warnings

- During the twelve weeks following surgery, the patient should avoid blows, vibrations, pressure and stress on the surgical site, and should not engage in any activity that could expose him/her to such risks.
- If any of the devices should break, it is recommended that the damaged devices be replaced if bone consolidation is not complete yet. Damaged components may migrate from their original location and/or injure the dura mater and brain. In any case, the final decision on whether to replace damaged devices lies with the surgeon.

Warning: You must report any serious incident that occurs in relation to the medical device supplied by us to the manufacturer, the Notified Body and the competent authority of the Member State in which you are established.

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